

Cruise ID: <u>2023-DIO-NOV</u>		Cast/Station Name: <u>all</u>		Analysis Date (dd mmm yyyy):	
Autosal Model: <u>8400B</u>		Serial Number: <u>73274</u>		<u>09 Nov 2023</u>	
Initial Standby Value: <u>24 + 5977</u>		Ratio After Standardization: <u>no change</u>			
Cell Standardization (Daily)	Standardization Values:	Standard Info.		Bath Temperature: <u>24</u> °C	
Conductivity must be ± 0.00001 compared to Standard	Before: _____	K ₁₅ VALUE: 0.999 <u>87</u>		Sample Temperature: <u>20.0</u> °C	
	After: _____	Batch #: <u>P166</u>		Room Temperature: <u>21.6</u> °C	
	Notes: _____	Batch Date: <u>exp. 6 April 2025</u>		Analyst: <u>ER</u>	

Notes: Run with -O80 samples.
These run first.
DNR not used in data processing - bottle was open.

